

Integrated OR at Elisabeth-TweeSteden - Tilburg, Netherlands

The Elisabeth-TweeSteden Hospital in Tilburg is one of the largest top clinical training hospitals in the Netherlands.

With a regional draw of almost half a million residents, in addition to patients from the North Brabant region, people from all over the Netherlands and neighboring Belgium come to this facility for medical care

Requirements

The Elisabeth-TweeSteden Hospital wanted a centralized workstation that could handle OR and image management according to the hospital's specific requirements.



"EIZO'S CuratOR solution offers maximum functionality and flexibility that we have completely integrated into our Operation Room."

Jan Steenbakker - Head Operating Room and Anesthesia



Solution

In close cooperation between the project managers and EIZO, a solution was developed that met all the hospital's requirements.

Customized for its specific OR, a wall console with four monitors was designed that not only visually fits in very well with the overall image, but also meets all functional requirements.

In the middle, both 24-inch monitors offer the option of PACS image viewing. Both 23-inch monitors outside allow the user to operate the OR management and control software – through the comfort of touchscreens. The connection to the PC and video management components to the external technology room is established by the customer-specific configuration of technology inside the Surgical Panel.

More information: www.eizo-or.com