



Operating Room Solutions

CuratOR™ Alipe

IP-based OR Infrastructure

EXPERIENCE THE FUTURE
OF OPERATING ROOM
TECHNOLOGY

EN





EVERYTHING THROUGH ONE CABLE – FORWARD-LOOKING OR INTEGRATION VIA IP

The future in the OR has already begun. Using EIZO's IP-based OR infrastructure, you can already intelligently link various devices and technologies. EIZO CuratOR Alipe expands the CuratOR Portfolio with a standardized, modular, expandable network solution that ensures the transmission of various data without discernible delay. Even over wide distances.

CuratOR Alipe is flexibly scalable – new components can simply be added as needed using plug and play.


A network-based infrastructure means infinite flexibility in the hospital and beyond – all data are always available everywhere.

Modalities in the OR
Connecting modalities, even through various interfaces



Mobile systems

The standardized network interface simplifies the integration of mobile systems



BENEFITS

- ✕ Universal transport medium: Transmission of video data, audio, control signals, USB, etc. via a data network
- ✕ Use of standardized technology: Backward compatibility, ongoing development, and expandability are guaranteed
- ✕ Designed for long distances
- ✕ Scalable and modular and therefore flexible with regard to future enhancements



Vast distances

Transmissions to distant rooms for training classes and conferences




TECHNICAL FACTS TRANSCIVER

- ✕ Encoder and decoder in one device, parameters set via software
- ✕ Uncompressed, loss-free data transmission via fiber-optic cable
- ✕ Video transmission of up to 10 GB / sec.
- ✕ Minimal delay for optimum hand-eye coordination
- ✕ Up to two video links can be processed simultaneously per device
- ✕ Multiple views possible on a monitor: Side-by-Side, Picture-in-Picture, QuadView
- ✕ Mouse and keyboard function via the network
- ✕ Audio and USB support
- ✕ Passive cooling; no resulting noise, no mechanical wear

◆ EIZO Locations

◆ EIZO Dealer Network



CONTACT PARTNER

You will find your EIZO contact partner at:
www.eizo-or.com/contact



EIZO GmbH
 Siemensallee 84
 76187 Karlsruhe, Germany

EIZO GmbH
 OR Solutions
 Branch Plauen
 Kopernikusstraße 80
 08523 Plauen, Germany

Tel: +49-3741-42365-0
 Fax: +49-3741-42365-199
 E-Mail: info@eizo-or.com
 Web: www.eizo-or.com

All product names are trademarks or registered trademarks of their respective companies.
 CuratOR is a trademark and EIZO is a registered trademark of EIZO Corporation.
 Specifications are subject to change without notice. Copyright © 2017 EIZO GmbH. All rights reserved.. (1704)